

IMPORTANT

This document is intended for healthcare professionals.

According to the Medical Device Directive 93/42/EEC, mammary implants are Class III and sizers for mammary implants are Class IIa medical devices. They are intended to be used in the context of plastic, reconstructive and aesthetic surgery. GROUPE SEBBIN reserves the use of its devices for physicians trained in the techniques of plastic, reconstructive and aesthetic surgery.

Please read the instructions for use carefully before use. They are available in private access on www.mysebbin.com.

The SEBBIN mammary implants and sizers for mammary implants are CE marked by the notified body mdc medical device certification GmbH number 0483.

Dimensions of our implants may present slight variations to the ones mentioned into this document.

*GROUPE SEBBIN follows the Standard FR EN ISO 14607: 2018 for the surfaces classification of its mammary implants.

MANUFACTURER



GROUPE SEBBIN SAS

39, Parc d'Activités des Quatre Vents
95650 Boissy-l'Aillier France
Tel.: +33 1 34 42 13 28
Fax: +33 1 34 42 16 88
www.sebbin.com

INTERNATIONAL OPERATIONS

Groupe Sebbin SAS, France

Tel.: +33 1 34 42 13 28 - Fax: +33 1 34 42 16 88
Email: serviceclients@sebbin.com

Sebbin Benelux SPRL

Tel.: +32 2 792 57 76 - Fax: +32 2 792 57 73
Email: info@sebbin.be

Sebbin Deutschland GmbH

Tel.: +49 2102 70 88 00-23 - Fax: +49 2102 70 88 00-99
Email for Germany: service@sebbin.de
Email for Austria: austria@sebbin.de

Sebbin Ibérica Productos Médicos, S.L.

Tel.: +34 93 159 3003 - Fax: +34 93 676 3053
Email: info@sebbin.es

Sebbin Suisse Sarl

Tel.: +41 43 336 72 30 - Fax: +41 43 336 72 36
Email: info@sebbin.ch

Sebbin UK Ltd

Tel.: 0800 652 4028 - Fax: +4433 0088 9575
Email: custserv@sebbin.uk

SALES CATALOGUE


VOLUPTY mammary implants

ROUND SHAPE MICRO-TEXTURED* SURFACE



Round mammary implants range

	PURITY SMOOTH*	SUBLIMITY SEMI-SMOOTH*	VOLUPTY MICRO-TEXTURED*
NATURGEL 1 - COHESIVE			
Soft consistency implants	LS 50 Moderate profile	LS 74 Semi-moderate profile	LS 90 Moderate profile
		LS 70 Moderate profile	
	LS 51 High profile	LS 71 High profile	LS 91 High profile
NATURGEL 2 - HIGH COHESIVE			
Firmer consistency implants	LSC 54 Moderate profile	LSC 76 Semi-moderate profile	LSC 92 Moderate profile
		LSC 72 Moderate profile	
	LSC 55 High profile	LSC 73 High profile	LSC 93 High profile



The arc is the length from the highest point of projection to the implant base.

VOLUPTY MICRO-TEXTURED* MAMMARY IMPLANTS
Naturgel 1 Cohesive

Moderate profile

Vol. (ml)	Base (mm)	Proj. (mm)	Arc (mm)	Reference	Sizer**
80	80	25	54	LS 90 080	GAB 35 080
100	84	28	59	LS 90 100	GAB 35 100
120	87	30	63	LS 90 120	GAB 35 120
140	91	31	66	LS 90 140	GAB 35 140
160	95	32	68	LS 90 160	GAB 35 160
180	101	33	72	LS 90 180	GAB 35 180
200	104	34	74	LS 90 200	GAB 35 200
220	108	35	75	LS 90 220	GAB 35 220
240	110	36	77	LS 90 240	GAB 35 240
260	112	38	78	LS 90 260	GAB 35 260
280	116	39	80	LS 90 280	GAB 35 280
300	119	40	82	LS 90 300	GAB 35 300
330	123	41	84	LS 90 330	GAB 35 330
360	126	42	87	LS 90 360	GAB 35 360
390	128	44	89	LS 90 390	GAB 35 390
430	133	45	93	LS 90 430	GAB 35 430
470	137	47	96	LS 90 470	GAB 35 470
510	140	49	98	LS 90 510	GAB 35 510
560	144	50	102	LS 90 560	GAB 35 560
610	147	51	108	LS 90 610	GAB 35 610



VOLUPTY MICRO-TEXTURED* MAMMARY IMPLANTS
Naturgel 1 Cohesive

High profile

Vol. (ml)	Base (mm)	Proj. (mm)	Arc (mm)	Reference	Sizer**
85	68	35	52	LS 91 085	GAB 36 085
105	72	39	55	LS 91 105	GAB 36 105
125	78	40	58	LS 91 125	GAB 36 125
145	82	41	62	LS 91 145	GAB 36 145
165	86	43	65	LS 91 165	GAB 36 165
185	91	44	67	LS 91 185	GAB 36 185
205	93	45	70	LS 91 205	GAB 36 205
225	97	46	73	LS 91 225	GAB 36 225
245	99	47	75	LS 91 245	GAB 36 245
265	101	49	78	LS 91 265	GAB 36 265
285	103	50	80	LS 91 285	GAB 36 285
305	105	52	82	LS 91 305	GAB 36 305
335	109	54	85	LS 91 335	GAB 36 335
365	112	56	87	LS 91 365	GAB 36 365
395	114	57	89	LS 91 395	GAB 36 395
435	117	60	91	LS 91 435	GAB 36 435
475	122	61	95	LS 91 475	GAB 36 475
515	126	62	99	LS 91 515	GAB 36 515
565	129	63	102	LS 91 565	GAB 36 565
615	134	65	105	LS 91 615	GAB 36 615



VOLUPTY MICRO-TEXTURED* MAMMARY IMPLANTS
Naturgel 2 High Cohesive

Moderate profile

Vol. (ml)	Base (mm)	Proj. (mm)	Arc (mm)	Reference	Sizer**
90	78	26	54	LSC 92 090	GAB 31 090
110	82	29	58	LSC 92 110	GAB 31 110
130	86	33	63	LSC 92 130	GAB 31 130
155	89	34	66	LSC 92 155	GAB 31 155
175	93	36	69	LSC 92 175	GAB 31 175
200	97	37	72	LSC 92 200	GAB 31 200
220	101	39	74	LSC 92 220	GAB 31 220
240	105	39	76	LSC 92 240	GAB 31 240
265	107	42	79	LSC 92 265	GAB 31 265
285	109	43	80	LSC 92 285	GAB 31 285
310	113	44	82	LSC 92 310	GAB 31 310
330	117	45	84	LSC 92 330	GAB 31 330
365	121	46	87	LSC 92 365	GAB 31 365
395	124	48	89	LSC 92 395	GAB 31 395
430	126	50	92	LSC 92 430	GAB 31 430
475	130	51	95	LSC 92 475	GAB 31 475
520	134	52	98	LSC 92 520	GAB 31 520
560	137	54	101	LSC 92 560	GAB 31 560
615	141	57	106	LSC 92 615	GAB 31 615
670	142	58	111	LSC 92 670	GAB 31 670



VOLUPTY MICRO-TEXTURED* MAMMARY IMPLANTS
Naturgel 2 High Cohesive

High profile

Vol. (ml)	Base (mm)	Proj. (mm)	Arc (mm)	Reference	Sizer**
90	66	39	55	LSC 93 090	GAB 32 090
110	70	41	58	LSC 93 110	GAB 32 110
130	75	42	60	LSC 93 130	GAB 32 130
150	80	43	64	LSC 93 150	GAB 32 150
175	84	45	67	LSC 93 175	GAB 32 175
195	87	48	71	LSC 93 195	GAB 32 195
215	91	50	73	LSC 93 215	GAB 32 215
235	94	51	75	LSC 93 235	GAB 32 235
255	96	53	78	LSC 93 255	GAB 32 255
280	98	55	80	LSC 93 280	GAB 32 280
300	101	56	82	LSC 93 300	GAB 32 300
320	103	57	84	LSC 93 320	GAB 32 320
350	106	58	86	LSC 93 350	GAB 32 350
385	109	60	89	LSC 93 385	GAB 32 385
415	112	62	91	LSC 93 415	GAB 32 415
455	116	64	94	LSC 93 455	GAB 32 455
495	119	66	99	LSC 93 495	GAB 32 495
540	122	67	102	LSC 93 540	GAB 32 540
590	127	68	105	LSC 93 590	GAB 32 590
645	132	69	108	LSC 93 645	GAB 32 645



**The sizers are prefilled with silicone gel with a smooth shell. They are delivered sterile, for single-use.