

MANUFACTURER



GROUPE SEBBIN SAS
39, Parc d'Activités des Quatre Vents
95650 Boissy-l'Aillie France
Tel.: +33 1 34 42 13 28
Fax: +33 1 34 42 16 88
www.sebbin.com

INTERNATIONAL OPERATIONS

Groupe Sebbin SAS, France

Tel.: +33 1 34 42 13 28 - Fax: +33 1 34 42 16 88
Email: serviceclients@sebbin.com

Sebbin Benelux SPRL

Tel.: +32 2 792 57 76 - Fax: +32 2 792 57 73
Email: info@sebbin.be

Sebbin Deutschland GmbH

Tel.: +49 2102 70 88 00-23 - Fax: +49 2102 70 88 00-99
Email for Germany: service@sebbin.de
Email for Austria: austria@sebbin.de

Sebbin Ibérica Productos Médicos, S.L.

Tel.: +34 93 159 3003 - Fax: +34 93 676 3053
Email: info@sebbin.es

Sebbin Suisse Sarl

Tel.: +41 43 336 72 30 - Fax: +41 43 336 72 36
Email: info@sebbin.ch

Sebbin UK Ltd

Tel.: 0800 652 4028 - Fax: +4433 0088 9575
Email: custserv@sebbin.uk

GS205-A-V04-2020/11 - Photo credit: Shutterstock - Creation: Sharing



IMPORTANT

This document is intended for healthcare professionals.

According to the Medical Device Directive 93/42/EEC, mammary implants are Class III and sizers for mammary implants are Class IIa medical devices. They are intended to be used in the context of plastic, reconstructive and aesthetic surgery. GROUPE SEBBIN reserves the use of its devices for physicians trained in the techniques of plastic, reconstructive and aesthetic surgery.

Please read the instructions for use carefully before use. They are available in private access on www.mysebbin.com.



The SEBBIN mammary implants and sizers for mammary implants are CE marked by the notified body mdc medical device certification GmbH number 0483.

Dimensions of our implants may present slight variations to the ones mentioned into this document.

*GROUPE SEBBIN follows the Standard FR EN ISO 14607: 2018 for the surfaces classification of its mammary implants.

Round mammary implants range

Table with 3 columns: PURITY SMOOTH*, SUBLIMITY SEMI-SMOOTH*, VOLUPTY MICRO-TEXTURED*. Rows include NATURGEL 1 - COHESIVE (Soft consistency implants: LS 50, LS 51, LS 74, LS 71, LS 90, LS 91) and NATURGEL 2 - HIGH COHESIVE (Firmer consistency implants: LSC 54, LSC 55, LSC 76, LSC 73, LSC 72, LSC 92, LSC 93).



The arc is the length from the highest point of projection to the implant base.



SALES CATALOGUE

SUBLIMITY mammary implants

ROUND SHAPE SEMI-SMOOTH* SURFACE



SUBLIMITY SEMI-SMOOTH* MAMMARY IMPLANTS
Naturgel 1 Cohesive

Semi-moderate profile

Vol. (ml)	Base (mm)	Proj. (mm)	Arc (mm)	Reference	Sizer**
80	82	21	56	LS 74 080	GAB 33 080
100	87	23	59	LS 74 100	GAB 33 100
120	92	24	62	LS 74 120	GAB 33 120
140	95	26	65	LS 74 140	GAB 33 140
160	103	27	69	LS 74 160	GAB 33 160
180	107	28	71	LS 74 180	GAB 33 180
200	111	29	73	LS 74 200	GAB 33 200
220	114	30	76	LS 74 220	GAB 33 220
245	118	31	78	LS 74 245	GAB 33 245
270	122	32	80	LS 74 270	GAB 33 270
295	126	33	83	LS 74 295	GAB 33 295
325	130	34	85	LS 74 325	GAB 33 325
370	136	35	88	LS 74 370	GAB 33 370
415	142	36	91	LS 74 415	GAB 33 415
470	148	38	94	LS 74 470	GAB 33 470
510	153	39	98	LS 74 510	GAB 33 510
560	157	42	102	LS 74 560	GAB 33 560
-	-	-	-	-	-
-	-	-	-	-	-
-	-	-	-	-	-



SUBLIMITY SEMI-SMOOTH* MAMMARY IMPLANTS
Naturgel 1 Cohesive

Moderate profile

Vol. (ml)	Base (mm)	Proj. (mm)	Arc (mm)	Reference	Sizer**
80	80	25	54	LS 70 080	GAB 35 080
100	84	28	59	LS 70 100	GAB 35 100
120	87	30	63	LS 70 120	GAB 35 120
140	91	31	66	LS 70 140	GAB 35 140
160	95	32	68	LS 70 160	GAB 35 160
180	101	33	72	LS 70 180	GAB 35 180
200	104	34	74	LS 70 200	GAB 35 200
220	108	35	75	LS 70 220	GAB 35 220
240	110	36	77	LS 70 240	GAB 35 240
260	112	38	78	LS 70 260	GAB 35 260
280	116	39	80	LS 70 280	GAB 35 280
300	119	40	82	LS 70 300	GAB 35 300
330	123	41	84	LS 70 330	GAB 35 330
360	126	42	87	LS 70 360	GAB 35 360
390	128	44	89	LS 70 390	GAB 35 390
430	133	45	93	LS 70 430	GAB 35 430
470	137	47	96	LS 70 470	GAB 35 470
510	140	49	98	LS 70 510	GAB 35 510
560	144	50	102	LS 70 560	GAB 35 560
610	147	51	108	LS 70 610	GAB 35 610



SUBLIMITY SEMI-SMOOTH* MAMMARY IMPLANTS
Naturgel 1 Cohesive

High profile

Vol. (ml)	Base (mm)	Proj. (mm)	Arc (mm)	Reference	Sizer**
85	68	35	52	LS 71 085	GAB 36 085
105	72	39	55	LS 71 105	GAB 36 105
125	78	40	58	LS 71 125	GAB 36 125
145	82	41	62	LS 71 145	GAB 36 145
165	86	43	65	LS 71 165	GAB 36 165
185	91	44	67	LS 71 185	GAB 36 185
205	93	45	70	LS 71 205	GAB 36 205
225	97	46	73	LS 71 225	GAB 36 225
245	99	47	75	LS 71 245	GAB 36 245
265	101	49	78	LS 71 265	GAB 36 265
285	103	50	80	LS 71 285	GAB 36 285
305	105	52	82	LS 71 305	GAB 36 305
335	109	54	85	LS 71 335	GAB 36 335
365	112	56	87	LS 71 365	GAB 36 365
395	114	57	89	LS 71 395	GAB 36 395
435	117	60	91	LS 71 435	GAB 36 435
475	122	61	95	LS 71 475	GAB 36 475
515	126	62	99	LS 71 515	GAB 36 515
565	129	63	102	LS 71 565	GAB 36 565
615	134	65	105	LS 71 615	GAB 36 615



SUBLIMITY SEMI-SMOOTH* MAMMARY IMPLANTS
Naturgel 2 High Cohesive

Semi-moderate profile

Vol. (ml)	Base (mm)	Proj. (mm)	Arc (mm)	Reference	Sizer**
80	82	21	56	LSC 76 080	GAB 30 080
100	87	23	59	LSC 76 100	GAB 30 100
120	92	24	62	LSC 76 120	GAB 30 120
140	95	26	65	LSC 76 140	GAB 30 140
160	103	27	69	LSC 76 160	GAB 30 160
180	107	28	71	LSC 76 180	GAB 30 180
200	111	29	73	LSC 76 200	GAB 30 200
220	114	30	76	LSC 76 220	GAB 30 220
245	118	31	78	LSC 76 245	GAB 30 245
270	122	32	80	LSC 76 270	GAB 30 270
295	126	33	83	LSC 76 295	GAB 30 295
325	130	34	85	LSC 76 325	GAB 30 325
370	136	35	88	LSC 76 370	GAB 30 370
415	142	36	91	LSC 76 415	GAB 30 415
470	148	38	94	LSC 76 470	GAB 30 470
510	153	39	98	LSC 76 510	GAB 30 510
560	157	42	102	LSC 76 560	GAB 30 560
-	-	-	-	-	-
-	-	-	-	-	-
-	-	-	-	-	-



SUBLIMITY SEMI-SMOOTH* MAMMARY IMPLANTS
Naturgel 2 High Cohesive

Moderate profile

Vol. (ml)	Base (mm)	Proj. (mm)	Arc (mm)	Reference	Sizer**
90	78	26	54	LSC 72 090	GAB 31 090
110	82	29	58	LSC 72 110	GAB 31 110
130	86	33	63	LSC 72 130	GAB 31 130
155	89	34	66	LSC 72 155	GAB 31 155
175	93	36	69	LSC 72 175	GAB 31 175
200	97	37	72	LSC 72 200	GAB 31 200
220	101	39	74	LSC 72 220	GAB 31 220
240	105	39	76	LSC 72 240	GAB 31 240
265	107	42	79	LSC 72 265	GAB 31 265
285	109	43	80	LSC 72 285	GAB 31 285
310	113	44	82	LSC 72 310	GAB 31 310
330	117	45	84	LSC 72 330	GAB 31 330
365	121	46	87	LSC 72 365	GAB 31 365
395	124	48	89	LSC 72 395	GAB 31 395
430	126	50	92	LSC 72 430	GAB 31 430
475	130	51	95	LSC 72 475	GAB 31 475
520	134	52	98	LSC 72 520	GAB 31 520
560	137	54	101	LSC 72 560	GAB 31 560
615	141	57	106	LSC 72 615	GAB 31 615
670	142	58	111	LSC 72 670	GAB 31 670



SUBLIMITY SEMI-SMOOTH* MAMMARY IMPLANTS
Naturgel 2 High Cohesive

High profile

Vol. (ml)	Base (mm)	Proj. (mm)	Arc (mm)	Reference	Sizer**
90	66	39	55	LSC 73 090	GAB 32 090
110	70	41	58	LSC 73 110	GAB 32 110
130	75	42	60	LSC 73 130	GAB 32 130
150	80	43	64	LSC 73 150	GAB 32 150
175	84	45	67	LSC 73 175	GAB 32 175
195	87	48	71	LSC 73 195	GAB 32 195
215	91	50	73	LSC 73 215	GAB 32 215
235	94	51	75	LSC 73 235	GAB 32 235
255	96	53	78	LSC 73 255	GAB 32 255
280	98	55	80	LSC 73 280	GAB 32 280
300	101	56	82	LSC 73 300	GAB 32 300
320	103	57	84	LSC 73 320	GAB 32 320
350	106	58	86	LSC 73 350	GAB 32 350
385	109	60	89	LSC 73 385	GAB 32 385
415	112	62	91	LSC 73 415	GAB 32 415
455	116	64	94	LSC 73 455	GAB 32 455
495	119	66	99	LSC 73 495	GAB 32 495
540	122	67	102	LSC 73 540	GAB 32 540
590	127	68	105	LSC 73 590	GAB 32 590
645	132	69	108	LSC 73 645	GAB 32 645



**The sizers are prefilled with silicone gel with a smooth shell. They are delivered sterile, for single-use.